MANTA

BODY FAT SCALE

MM442P

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3. Selection buttons

4. ball-of-foot electrodes

PREPARATION BEFORE USE*Inserting the Batteries**

Open the battery cover on the back of the measuring platform. Pay attention to the Anode and Cathode, then put the Battery into the battery Compartment.

Note: If you do not intend to use this unit for a long period of time, it is advisable to remove the battery before storage.

POSITIONING THE MONITOR

Place the measuring platform on a hard, flat surface where there is minimal vibration to ensure safe and accurate measurement.

Note: To avoid possible injury, do not step on the edge of the platform.

***HANDLING TIPS**

This monitor is a precision instrument utilizing state of the art technology. To keep the unit in the best condition, follow these instructions carefully.

- Do not attempt to disassemble the measuring platform.
- Place the unit in an area free from direct sunlight, heating equipment, high humidity, or extreme temperature change.
- Never submerge in water. Use alcohol to clean the electrodes and glass cleaner(applied to cloth first) to keep them shiny, soaps is not allowed.
- Do not step on the platform when wet.
- People with pacemakers can not use weights.
- People with metal implants fail to obtain accurate results (not results of measurement of body weight).
- In order to ensure comparability of measurements we recommend that abywykonywać measurements always at the same time.

***USEFUL FEATURES**

- This Body fat Scale allow you to measure weight, body fat percentage, water percentage, muscle percentage and bone weight simultaneously and easily simply by stepping on the scale.
- Personal data(Female/Male, Age, Height) can be pre-set and stored in personal data memories.
- After you have completed the measuring process, the power shuts off automatically .
- With Automatic Zero function.
- With Low battery and overload indication function.

***WEIGHT ONLY FUNTION**

This scale is with Weight Only function (auto on) when the user doesn't press any key to set any data , only step on the plafrom, the display will show the user as weight only (the weight will flash twice and hold finally).

*HOW TO DETERMINE BODY FAT PERCENTAGE

Setting and Storing Data in Memory The unit can be operated if data has been programmed into one of the personal data memories

1. Turn on the Power

Press SET button to turn on the unit. The unit will confirm activation; the last user's personal data will be appeared alternately.

2. Use the Selector Switch on the bottom to select a weight scale to use: Metric(kg), Pounds(lb), or Stones(ST).

3. Select a Personal data Code

Firstly Press the SET button, the Personal data Code will be flashed. Press the Up/Down buttons to select a Personal data code. Once you reach the Personal data code you wish to use, press the set button to confirm.

4. Select Female or Male

Use the button Up/Down buttons to scroll through female, male settings, and then press the Set button to confirm.

5. Specify the Height

Use the Up/Down buttons to specify Height and then press the Set button to confirm.(range of user height is between 100cm-220cm)

6. Set Age

Use the Up/Down buttons to specify age. When you reach your age, press the Set button to confirm.(range of user age is between 10-80 year)

NOTE: If you make a mistake before you have finished programming it, repeat the steps from Select a Personal data Code.

*TAKING WEIGHT AND BODY FAT READINGS

After programming your personal data, you are ready to take a reading. Select Personal Data Code Press the Up button turn on. Use the Up/Down buttons to select your Personal data code. The display will show the programmed data. Until the display will be show "0.0", now step onto the platform.

Note: To ensure accuracy reading should be with barefoot to step on the platform and be sure that your heels are correctly aligned with the electrodes on the measuring platform. If you step onto the platform before "0.0" appears, you will not obtain a reading. Furthermore, if you do not step onto the measuring platform within about 10 seconds after "0.0" appears, the power is shut-off automatically.

Get Your Reading

Your weight will be shown first. Continue to stand on the platform "0000" will appear on the display and disappear one by one from left to right. The display will then flash your weight, Then the body fat percentage water, percentage bone weight, muscle percentage, calorie and weight will appear alternately twice.

Note:

Do not step off from until when body fat is shown. If the body fat is shown, even if you step off the platform. The alternating display, body fat percentage and water percentage, bone weight, muscle percentage, calories and weight will appear alternately twice, the unit will then shut-off automatically.

Body Fat Ratio Chart (%)

Age	FEMALE (%)				MALE (%)			
	Low	Health	High	Obese	Low	Health		
10-39	14-20	21-25	26-31	32-38	11-17	18-23	24-29	30-36
40-55	15-21	22-26	27-32	33-38	12-18	19-24	25-30	31-37
56-80	16-22	23-27	28-33	34-38	13-19	20-25	26-31	32-38

Water ratio chart (%)

Age	FEMALE (%)			MALE (%)		
	Low	Healthy	High	Low	Healthy	High
10-15	<57	57-67	>67	<58	58-72	>72
16-30	<47	47-57	>57	<53	53-67	>67
31-60	<42	42-52	>52	<47	47-61	>61
61-80	<37	37-47	>47	<42	42-56	>56

Weight conversion table

1 lb (pound)	=	0,45 kg (kilogram)
1 kg (kilogram)	=	2,2 lb (pound)
1 st (Stone)	=	6,35 (kilogram)

The demand for calories:

Note!

The data displayed is an estimate of the minimum value needed for the body to maintain specified when measuring the user's weight.

* Establishment of the passive live mode

The real value of total calories should be selected individually, it depends on factors such as physical activity and metabolism. For professional advice you should go to a doctor or nutritionist.

***REPLACING THE BATTERIES**

When the batteries run low, the Lo message will appear on the display window. When the message appears, be sure to replace the batteries immediately. Since weak batteries will affect the accuracy of your measurements. Change all the batteries at the same time with new batteries.

Note: Your settings will not be erased from the memory when you remove the Batteries.

***ERROR MESSAGE**

The following error messages may appear:

-Lo: The battery in the scales must be replaced.

-Err: The maximum load capacity is exceeded.(range of weight with 3-151kg body fat fuction:0.8-150kg))

-Err2%: The body-fat and body-water percentages cannot be calculated. (If body fat percentage is less than 5% or more than 85%, the reading cannot be accurate.)

Note!

Product design and Technical parameters may be change without notice. This mainly concerns parameters technical, software and user manual present. User's Manual is a general orientation on service product. Manufacturer and distributor do not assume any responsibility as compensation for any inaccuracies errors in the descriptions appearing in this user manual.

Specification

Voltage:	20-240V~50/60Hz
Touch Buttons:	Tak
Battery included:	2x1,5V AAA
Scale bar:	100 g
Accuracy:	100 g
Units of Measurement:	Kilograms (kg), Pounds (lb) Stones (ST).
Maximum weight:	150kg
Additional measurements:	The percentage of fat, The percentage of water content, The content of calcium in the bones Percentage of muscle tissue Daily requirement of calories Weight.
Calibration:	Automatic:
Number of Pre users	10
Type of weight:	Electronic
Reading of measurement:	LCD screen
Platform:	Flat

INFORMATION ON THE CORRECT PROCEEDINGS IN THE CASE OF DISPOSING OF ELECTRONIC AND ELECTRIC DEVICE

The crossed dust bin symbol placed on our products indicates that the redundant electric or electronic products should not be disposed together with household waste.

Appropriate action in the event of having to dispose of equipment or components, or recycle them involves bringing the machine to a collection point where it will be accepted free of charge. In some countries you may dispose of the product when purchasing a new device.

The correct disposal of this product makes it possible to preserve the earth's natural resources for longer and prevents the degradation of the environment. Information about collection points is issued by local authorities. Incorrect disposal of waste is subject to penalties provided for in the law of the relevant area.

CUSTOMER SERVICE

tel: +48 22 332 34 63 or e-mail: serwis@manta.com.pl

FROM MONDAY TO FRIDAY, IN HOURS: 9.00-17.00

The manufacturer reserves the right to introduce changes