# MANTA

## ELECTRONIC BABY SCALE

MM443B

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## Spis treści

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Made in P.R.C.  
**FOR MANTA EUROPE**

Dear customer,

Thank you for choosing a product of our range. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the information they contain.

### **CARE AND ATTENTION**

- In normal weighing, the scale only receives static pressure. If the scale is used abnormally by impact, or dropped pressure, the glass may break easily and may cause danger.
- Please stand evenly on the scale and do not move.
- Do not use any chemical cleaning agent. Clean your scale with damp cloth, but do not allow water to get inside the scale.
- Put the scale in a cool and dry place. Do not store it in wet environment. As it may destroy the battery and shorten the battery life.
- Put the scale on flat surface; don't put it in an upright position.
- Remove the battery if the scale is not going to be used for a long period of time or in transportation.
- If there is no display after the scale starts, please check whether the battery and battery contacts are well connected. If the scale has not been used for a long period, please check that you need to replace the battery or not.
- The scale is for family weighing use only not for commercial or other purposes.
- Please do not put heavy objects on the scale.

### **CAUTION**

- Protect your personal scales from impact with hard objects, moisture, dust, chemicals, toiletries, liquid cosmetics, great temperature fluctuation and closeness to sources of heat (open fires, radiators).
- Cleaning: You can clean the scale with a damp cloth and a little washing up liquid, if required. Never immerse the scale in water or rinse it under running water.
- Storage: Do not place any objects on the scale when not in use.
- Strong electromagnetic fields (e.g. cell phones) may impair the accuracy of the scale.
- The scales can accept a maximum load of 44 lb (20 kg/704 oz). The weight measurement is displayed in 0.20 oz (5 g/0.02 lb) divisions.
- Remove transport safety catch.
- Place the scale on a level and solid surface; this is important for ensuring that the scale weighs correctly.

### **OPERATION**

#### **Battery**

The scale requires 1XCR2032 battery. Before using the baby scales, these must be inserted in the battery compartment beneath the scales, ensuring the right polarity. When "Lo" appears in the display, the battery need to be replaced.

**Weighing**

The weight is displayed in steps of 5 g from a minimum weight of 50 g.

Press the button "" to switch on the scale. As soon as "0.000 kg" appears on the display, the scale is ready to take measurements. Now place the baby on the weighing area. Do this very carefully; do not let the baby fall and take particular care with the baby's head.

This baby scale has a "HOLD" function (AUTO HOLD and MECHANICAL HOLD). The AUTO HOLD function establishes the baby's weight on the display as soon as he/she is still enough to achieve an exact measurement value. If the baby is moving around vigorously, his/her weight cannot be recorded. Therefore, try to calm the wriggling baby to achieve a rapid measurement, but do not touch the baby, as this will affect the weight measurement. As soon as the baby's weight is established on the display, the word "HOLD" stops flashing and the weight is displayed for a further 15 seconds. The scale then switches off automatically. If the baby is wriggling about too vigorously, you can select the MECHANICAL HOLD function by pressing the "HOLD" button. Within 3 seconds, the word "HOLD" stops flashing and the weight is displayed for a further 15 seconds. If you do not want to put the baby directly on the surface of the weighing, but prefer, for example, use a towel, you have to put a towel on weight before turning it on. Weight turn on automatically reset. TARE function can also be activated manually. To do this, turn on the weight and then put on weight all the extra items. To reset the scale to 0,000 kg, press the "TARE" when "HOLD" will stop flashing. The display shows the additional inscription TARA. At this point you can gently put the baby on the scale and start weighing process.

On delivery the scale is set at "kg" units. The button "kg/lb/oz" allows you to switch at any time between kilograms units and pounds (lb) and ounces (oz).

**SPECIFICATIONS**

- Equipped with a High Precision "Strain-Gauge" Sensor.
- Large LCD Display
- Tare weight function: Auto and Manual
- Max capacity:20kg
- Graduation:5g
- Adjustable KG/LB/OZ
- Automatic switch off
- Battery included

**Note!**

Product design and Technical parameters may be change without notice. This mainly concerns parameters technical, software and user manual present. User's Manual is a general orientation on service product. Manufacturer and distributor do not assume any responsibility as compensation for any inaccuracies errors in the descriptions appearing in this user manual.

**Note about disposal**

Older appliances, which display the symbol as shown in the picture, may not be discarded in household waste. You must take them to a specific disposal site (check with your local authority) or to your retailer. These sites will guarantee disposal without harming the environment

### **CUSTOMER SERVICE**

tel: +48 22 332 34 63 lub e-mail: [serwis@manta.com.pl](mailto:serwis@manta.com.pl)  
FROM MONDAY TO FRIDAY, IN HOURS: 9.00-17.00

The manufacturer reserves the right to make changes to product specifications without notice.

Main page: [www.manta.com.pl](http://www.manta.com.pl)  
Technical suport page: [www.manta.info.pl](http://www.manta.info.pl)

Made in P.R.C.  
**FOR MANTA EUROPE**