MANTA

Thermometer

THM501P

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SPIS TREŚCI

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Made in P.R.C.
FOR MAN TA EUROPE

1. General Description

Non-contact Forehead IR Thermometer is specially designed to take the body temperature of a person (or Animal) regardless of room temperature. Depending on various skin types and thickness, there may be temperature difference.

2. Safety Information

- ✱ This device must only be used for the purposes described in this instruction manual..
- ✱ This device must only be used in an ambient temperature range between 10 and 40°C.
- ✱ Do not rely on temperature measurements alone, If you have health concerns, seek medical advice.
- ✱ Do not expose this thermometer to electric shocks.
- ✱ Do not expose this thermometer to extreme temperature conditions of > 50°C or < 0°C.
- ✱ Do not use the device in relative humidity higher than 85%.
- ✱ Do not use the device near large electromagnetic fields such as found with cordless or cell phones.
- ✱ Keep the device away from water and heat, including direct sunlight.
- ✱ Do not drop or knock the device and do not use if damaged.
- ✱ It may affect the accuracy of measurements when the forehead is covered by hair, perspiration, cap or scarf.(see part 10-4).
- ✱ Keep the Measuring distance as 5-15cm(2in-5.9in) , See part 10-4).
- ✱ When the body infrared thermometer should be left in that room during 15 to 20 minutes before using.
- ✱ It may affect the accuracy of measurements when the forehead is covered by perspiration or other factors, please take the temperature behind the ear lobe(See part 10-5).
- ✱ Clean the glass with a cotton bud lightly moistened with 70% alcohol.
- ✱ Remove the batteries if not in use for a long period.

Importance

- ✱ Before taking the temperature make sure to remove hair and perspiration from the forehead.
- ✱ Selecting "Body" mode to measure the body temperature, selecting "Surface" to measure the surface temperature.
(Selecting "Animal" mode to measure a animal, such as a cat or dog)
- ✱ Should a problem occur with your device, please contact your retailer. Do not attempt to repair it by yourself.
- ✱ According to EMC standard, the medical electronic products should be maintained specially
- ✱ Application: This proeluet is suitable for any peoples body temperature in any age.

3. Features

- ✱ Precise non-contact measurements
- ✱ User selectable °C or °F
- ✱ Set Alarm value
- ✱ Memorization of the last 99 measurements
- ✱ Automatic Data Hold and auto power off
- ✱ Automatic selection range and display resolution 0.1°C(0.1°F)
- ✱ 3 colors backlight LCD display
- ✱ 3 modes: Body mode, Surface mode , Animal mode.

4. Intended Use

Non-contact Forehead IR thermometer is designed for body surface and forehead temperature measurement for infants and adults without contact to human body.

It can also be used to measure the temperature of baby-bottle or both, or room temperature(by using the surface mode) and animal.

Normal Temperatures According to Measurement Method

Measurement Method	Normal Temp. °C	Normal Temp. °F
Rectal	36.6 to 38.0	97.8 to 100.4
Oral	35.5 to 37.5	95.9 to 99.5
Axillary	34.7 to 37.3	94.4 to 99.1
Ear	35.8 to 38.0	96.4 to 100.4

The temperature of the human body varies throughout the day. It can also be influenced by numerous external factors: age, sex, type and thickness of skin...

Normal Temperatures According to Age

Age	Temp. °C	Temp. °F
0-2 Years	36.4 to 38.0	97.5 to 100.4
3-10 Years	36.1 to 37.8	97.0 to 100.0
11-65 Years	35.9 to 37.6	96.6 to 99.7
> 65 Years	35.8 to 37.5	96.4 to 99.5

Normal Temperature for Animals


Animal	Temp. °C	Temp. °F
Cat	38.1 to 39.2	100.6 to 102.6
Dog	37.9 to 39.9	100.2 to 103.8
Rabbit	38.6 to 40.1	101.8 to 104.2
Pig	38.7 to 39.8	101.7 to 103.6
Sheep	38.3 to 40.2	100.9 to 104.4
Cow	36.7 to 39.5	98.1 to 103.1
Chook(Bird)	39.6 to 43.6	103.3 to 110.5
Goose	40.0 to 44.0	104.0 to 111.2
Duck	40.0 to 41.0	104.0 to 105.8

5. Configuration

1. IR sensor
2. LCD display
3. Mode selection
4. Down button
5. Up button
6. Mode button
7. Measurement Trigger
8. Battery cover

6. Indicator

1. Surface mode symbol
2. Body mode symbol
3. Animal mode
4. Digital readout
5. Battery symbol
6. Memory symbol
7. Save data readout
8. Temp. °C(Celsius)/
°F(Fahrenheit) scale
9. Buzzer symbol
10. Laser aiming

7. Description of Symbols

CE1023	The device is in accordance with Medical Device Directive 93/42/EEC
	The device is in accordance with FCC part 15 subpart B:2007/Radio Frequency Devices.
	3V DC power supply
	Type B equipment
	In order to protect the environment, please recycle the battery according to the local regulations.
Body /Surface /Animal	Indication of Mode
	Attention , consult accompanying documents
	This manual is made of paper, which can be used circularly.

8. Technical Specifications

Product Name	Multi-Purpose Non Contact Infrared Thermometer
Model	THM501P
Display Resolution	0,1°C(0,1°F)
Operating Temp.	10 to 40°C (50 to 104°F)
Storage Temp.	0 to 50°C (32 to 122°F)
Humidity Rate	≤85%
Power	3.0V (2X "AAA" batteries)
Size	149(L)×95(W)×45(D) MM/ 5.9×3.7×1.8 inch
Weight	140g
Measurement Range	Body 32.0 to 42.9°C(89.6 to 109.3°F)
	Surface 0 to 100°C(32 to 212.0°F)
	Animal 32.0 to 45.0°C(89.6 to 112.9°F)
Precision	Human ±0.2°C(0.4°F) 36.0 to 39.0°C(96.8 to 102.2°F)
	±0.3°C(0.5°F) 32.0 to 35.9°C (89.6 to 96.2°F)
	±0.3°C(0.5°F) 39.0 to 42.9°C(102.2 to 109.2°F)
	Surface ±1.0°C(1.8°F) 0 to 100°C(32 to 212.0°F)
	Animal ±0.4°C (0.7°F) 32.0 to 45.0°C(89.6 to 112.9°F)
Measuring Distance	5 to 15 CM/ 2 to 5.9 inch
Backlight	Body mode 32.0 to 37.4°C(86.0 to 99.3°F)Green 37.5 to 37.9°C(99.5 to 100.2°F)Orange 38.0 to 42.9°C(100.4 to 109.2°F)Red Surface mode and Animal mode always keep Green color
Automatic stop	7 Sec

9. Advice for The First Use

For stable and reliable results, it is essential to check the Non-contact infrared thermometer and change as needed, as follows:

1st step:	Take the temperature of a person using a conventional thermometer, you will get 37.5°C(99.5°F) for instance
2nd step:	Take the temperature of the same person using the device keeping the 5 to 15 cm (2 to 5.9 in) distance between the thermometer and the forehead(Take care to remove any obstacle which could alter the measurement(hair, perspiration...)). If you get 37.5°C (99.5°F), the device is properly set and ready for use. If you get a lower temperature, such as 36.4°C (97.5°F), your difference is 1.1°C(2.2°F). You should adjust the temperature on the device and add the difference, i.e.1.1°C(2.2°F). To do it, press the MODE button for 2 seconds, the screen displays F1, press MODE button again until you get F3, Press UP button in order to add the difference(in our example, 1.1°C(2.2°F).
3rd step:	To check, take the temperature again using the product.


10. Use

10-1 Install batteries

10-2 For the first use or when inserting new batteries, please wait about 10 minutes for the warm-up of the apparatus.

10-3 If the device is not used for a long time, once you turn it on again, the device will test the room temperature first and will delay turning on for one or two seconds.

10-4 Select body mode, aim towards the forehead(see the diagram below for the positioning), from a distance of 5cm(2in), press the measuring key, the temperature is displayed immediately, accompanied by a beep sound. Making sure there is no hair, perspiration, cosmetic or cap covered on the forehead.

10-5 When the room temperature is significant difference, or there is perspiration on the forehead, you can take the temperature behind the ear lobe, Making sure there is no hair, perspiration, cosmetic or cap covered.

If the measurement is 38.0°C or more(can be setting, refer to 11-2), the sound should be be-be-be-be-be 5 rapid sounds.

10-6 Select object mode for measuring room temperature, liquid temperature, ect.

10-7 Select animal mode for taking temperature of an animal.

11. Setting Operation

	Mode Code	Name	Operation
11-1	F1	Set °C/°F scale	Press MODE button for 2 seconds, the screen will display F1 on the bottom right, press UP or DOWN button can change or °C
11-2	F2	Set Alarm temp.	Press MODE button for 2 seconds, the screen display F1, press MODE button to get F2, Select UP to increase the threshold by 0.1°C(0.1°F), DOWN to reduce it by. The alarm threshold default value is 38.0°C(100.4°F)
11-3	F3	Set difference	Press MODE button for 2 seconds, enter into F1, Press two times to get F3, select UP to increase the difference by 0.1°C(0.1°F),DOWN to reduce it by 0.1°C(0.1°F), In the case of seasonal or environmental changes a verification and adjustment should be carried out. Note This function is only effective body.
11-4	F4	Set Buzzer on/off	Press MODE button for 2 seconds, the screen displays F1, Press 3 times MODE button to get F4. Select UP to turn on the buzzer, press DOWN to stop it.
11-5	F5	Set laser aiming on/off	Press MODE button for 2 seconds, the screen displays F1, Press 4 times MODE, enter into F5, Select UP or DOWN to open or close the laser aiming function.

11-6	Measurement Mode	There are 3 measurement modes for this device. That is Body, Object and Animal. It is able to take the body temperature of human beings, take the temperature of an area or an object, a food, a liquid or a room, take the temperature of an animal, such as a pig, a dog, a cat, a cow...using different mode.	
		Body mode:	Put the mode button on the left, Measurement range 32.0 to 42.9°C (86.0 to 109.2°F).
		Object mode:	Put the mode button in the middle, Measurement rang 0 to 100°C(32.0 to 212.0°F)
		Animal mode:	Put the mode button on the right, Measurement range 32.0 to 45.0°C(86.0 to 113.0°F)
		Note: If your device has only 2 measurement modes, it is normal. Please ignore animal mode and consult the saler.	
11-7	3 Colors Backlight	In body mode, the device is able to judge whether a fever or not by different color of backlight: 32.0 to 37.4°C (86.0 to 99.3°F) Green 37.5 to 37.9°C (99.5 to 100.2°F) Orange 38.0 to 42.9°C (100.4 to 109.2°F) Red In Object mode or Animal mode, the backlight always keeps green.	
11-8	Data Memory	Data memory automatically after temperature measurements, which will display at the right corner of LCD, Press UP or DOWN button to display the last temperature measurement. Under the condition of power off, press UP and DOWN button together for two seconds, to display the last temperature measurement. At this mode if change the order number to 0 and press and press MODE button, then delete all memory data.	

11-9	Replace the Batteries	When the LCD displays: "the battery is used up." Open the lid and replace the batteries, taking great care with the correct positioning, A mistake with this could cause damage to the apparatus and compromise the guarantee of your device. Never use rechargeable batteries. Use only batteries for single usage. Remove the battery from the instrument if it is not required for extended periods of time in order to avoid damage to the thermometer resulting from a leaking battery.
11-10	Longevity Use	The device was conceived for an intense and professional use, its longevity is guaranteed for 18 months.

12. Advice

- * The protective glass over the lens is the most important and fragile part of the thermometer, please take great care of it.
- * Do not recharge non rechargeable batteries, do not throw in fire.
- * Do not expose the thermometer to sunlight or water.

13. Care and Cleaning

- * The infrared sensor is the most delicate part of the device,
- * It has to be clean and intact to ensure accurate readings.
- * Clean the device gently with a cotton swab or soft cloth moistened with 70% alcohol
- * Use a soft, dry cloth to clean the thermometer display and exterior.
- * Do not use abrasive cleaners, Never submerge the thermometer in water or any other liquid.
- * Store the device in a dry environment, and keep it away from dust and direct sunlight.
- * Always keep the thermometer out of reach of children.

14. Trouble shooting

Trouble	Reason	Solution
No response	Batteries used up?	Replace new batteries
	Batteries in wrong polarity?	Adjust batteries to proper polarity
Lo, Hi	Exceed measurement range Body mode: $< 32.0^{\circ}\text{C}$ display Lo $> 42.9^{\circ}\text{C}$ display Hi Object mode: $< 0^{\circ}\text{C}$ display Lo $> 100^{\circ}\text{C}$ display Hi Animal mode $< 32.0^{\circ}\text{C}$ display Lo $> 45.0^{\circ}\text{C}$ display Hi	1. Make sure the probe is clean and there is no obstruction prior to taking a temperature. 2. Make sure there is no air flux as this could interfere with the infrared system. 3. Make sure the measuring distance is not too far.
Err	Ambient temperature too low or too high	Allow the device to remain in a room for 20 minutes where the temperature is range from 10 to 40°C
	The batteries are running low	Replace 2 new AAA batteries

15. This appliance conforms to the following standards

ASTM E1965-1998

- EN980: Graphical symbols for use in the labling of medical devices
- En1041: Information supplied by the manufacturer with medical devices
- EN60601-1: Medical electrical equipment part 1:General requirements for safety(IEC60601-1:1998)
- EN60601-1-2: Medical electrical equipment part 1-2: General requirements for safety collateral standard electromagnetic compatibility requirements and test(IEC60601-1-2:2001)

EMC statement

This device has been tested and homologated in accordance with EN60601-1-2 for EMC, This does not guarantee in any way that the device will not be affected by electromagnetic interference. Avoid using the device in high electromagnetic environment.

The medical delectrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the accompanying documents.

INFORMATION ON THE CORRECT PROCEEDINGS IN THE CASE OF DISPOSING OF ELECTRONIC AND ELECTRIC DEVICE

The crossed dust bin symbol placed on our products indicates that the redundant electric or electronic products should not be disposed together with household waste.

Appropriate action in the event of having to dispose of equipment or components, or recycle them involves bringing the machine to a collection point where it will be accepted free of charge. In some countries you may dispose of the product when purchasing a new device.

The correct disposal of this product makes it possible to preserve the earth's natural resources for longer and prevents the degradation of the environment. Information about collection points is issued by local authorities. Incorrect disposal of waste is subject to penalties provided for in the law of the relevant area.

**CUSTOMER SERVICE**

tel: +48 22 332 34 63 lub e-mail: serwis@manta.com.pl
FROM MONDAY TO FRIDAY, IN HOURS: 9.00-17.00

The manufacturer reserves the right to make changes to product specifications without notice.

Main page: www.manta.com.pl
Technical suport page: www.manta.info.pl

Made in P.R.C.
FOR MANTA EUROPE